
Product Data Sheet

Product Name: Boc-Ala-ONp

Cat. No.: GA20892

Chemical Properties

Cas. No. 2483-49-0

Formula $C_{14}H_{18}N_2O_6$ M.Wt 310.31

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for cathepsin B and porcine pancreatic and human leukocyte elastases. It is also hydrolyzed by chymotrypsin and trypsin as well as by the respective zymogens, especially chymotrypsinogen. The catalytic activity of both chymotrypsinogen and tr

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA