
Product Data Sheet

Product Name: Boc-AOAc-OH

Cat. No.: GA20896

Chemical Properties

Cas. No. 42989-85-5

Formula C₇H₁₃NO₅

M.Wt 191.18

Solubility DMSO : 100 mg/mL (523.07 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Protected aminooxyacetic acid derivative for the incorporation of an N-terminal O-alkylhydroxylamine moiety, which is used in the chemoselective ligation of unprotected peptide fragments, for the conjugation to lipids and glycosides, and for the introd

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA