
Product Data Sheet

Product Name: Boc-Homoarg(Et)₂-OH (symmetrical)

Cat. No.: GA21030

Chemical Properties

Cas. No. 122532-94-9

Formula C₁₆H₃₂N₄O₄ M.Wt 344.46

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

N ω , ω' -alkylated derivative of homoarginine. Building block for LHRH antagonists as Ganirelix, which contains both enantiomers of Har(Et).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA