
Product Data Sheet

Product Name: Bz-Arg-pNA . HCl

Cat. No.: GA21201

Chemical Properties

Cas. No. 21653-40-7

Formula $C_{19}H_{22}N_6O_4 \cdot HCl$

M.Wt 434.88

Solubility DMSO : 100 mg/mL (229.95 mM; Need ultrasonic); H₂O : 5 mg/mL (11.50 mM; ultrasonic and warming and heat to 60°C) Store Storage at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for trypsin. BApNA is also cleaved by papain, acrosin, and plasmin.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA