
Product Data Sheet

Product Name: CD36 (93-110)-Cys

Cat. No.: GA21252

Chemical Properties

Cas. No. 143257-75-4

Formula C₉₆H₁₅₁N₂₉O₃₃S

M.Wt 2271.5

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

The CD36 sequence YRVRFLAKENVVTQDAEDNC blocked binding of CD36 to immobilized thrombospondin and partially inhibited collagen-induced platelet aggregation.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA