
Product Data Sheet

Product Name: Checkpoint Protein Hus1 (213-232) (human)

Cat. No.: GA21261

Chemical Properties

Cas. No.

Formula C₉₅H₁₅₀N₃₀O₃₅S

M.Wt

2304.48

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

HDAC1 plays an important role in modulating the eukaryotic chromatin structure and was found to interact specifically in yeast, mammalian cells, and in vitro with the human Hus1 gene product, as reported by Cai and coworkers. The Hus1 protein of *S. pom*

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA