
Product Data Sheet

Product Name: Collagen Type I α 1 chain (502-516)

Cat. No.: GA21275

Chemical Properties

Cas. No.

Formula C₆₃H₉₅N₁₉O₂₁

M.Wt

1454.54

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

GF-Hyp-GERGVEGP-Hyp-GPA corresponds to an integrin α 2 β 1-binding site of collagen I. The GER motif is essential for recognition.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA