
Product Data Sheet

Product Name: Cys-CD36 (139-155)

Cat. No.: GA21350

Chemical Properties

Cas. No. 143257-74-3

Formula C₈₇H₁₃₃N₂₅O₂₆S

M.Wt

1977.23

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The CD36 glycoprotein IV or IIIb sequence CNLAVAAASHIYQNQFVQ represents part of the epitope that binds the monoclonal anti-CD36 anti-body OKM5. It enhanced binding of CD36 to thrombospondin and augmented ADP-induced and collagen-induced aggregation in

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA