
Product Data Sheet

Product Name: Decanoyl-Arg-Arg-Leu-Leu-chloromethylketone

Cat. No.: GA21370

Chemical Properties

Cas. No. 911379-57-2

Formula C₃₅H₆₇ClN₁₀O₅ M.Wt 743.43

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

A cell-permeable irreversible inhibitor of proprotein convertase site 1 protease (S1P, also called subtilisin kexin isozyme, SKI-1). Dec-RRL-CKM efficiently blocked the processing of Lassa virus glycoprotein precursor at nanomolar concentrations, thou

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA