
Product Data Sheet

Product Name: Decanoyl-Arg-Val-Arg-Lys-chloromethylketone

Cat. No.: GA21371

Chemical Properties

Cas. No. 534615-50-4

Formula C₃₄H₆₆ClN₁₁O₅

M.Wt

744.42

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

The inhibitor Dec-RVRK-CMK has been used to characterize the specificity of furin-like proteases.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA