
Product Data Sheet

Product Name: Fibrinogen γ -Chain (117-133)

Cat. No.: GA21456

Chemical Properties

Cas. No. 160927-63-9

Formula $C_{84}H_{147}N_{25}O_{27}$ M.Wt 1939.24Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This fibrinogen sequence binds to ICAM 1 and mediates leukocyte-endothelium interaction. It potentially contributes to fibrinogen dependent inflammation. This provides a target for highly selective inhibition of fibrinogen dependent cellular inflammation

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA