
Product Data Sheet

Product Name: Fibronectin CS-1 Fragment (1978-1982)

Cat. No.: GA21457

Chemical Properties

Cas. No. 150525-67-0

Formula $C_{26}H_{45}N_5O_{10}$ M.Wt 587.67

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

EILDV and EILDVPST are derived from the type III connecting segment domain (III_{CS}) of fibronectin and have been found to promote tumor cell attachment after surface immobilization and to inhibit the experimental tumor metastasis of B16-BL6 melanoma. He

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA