
Product Data Sheet

Product Name: Fibronectin Fragment (1376-1380)

Cat. No.: GA21458

Chemical Properties

Cas. No. 158622-13-0

Formula C₂₄H₃₉N₁₁O₈ M.Wt 609.64

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The short peptide PHSRN has been identified in the 9th type III fibronectin repeat as a key site for synergistic enhancement of the cell-adhesive activity of the RGD sequence in the 10th type III repeat. However, this synthetic peptide mediated a syner

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA