
Product Data Sheet

Product Name: FITC-εAhx-Antennapedia Homeobox (43-58) amide

Cat. No.: GA21463

Chemical Properties

Cas. No. 1118750-35-8

Formula $C_{131}H_{191}N_{37}O_{25}S_2$

M.Wt 2748.32

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

FITC-LC-Antp Homeobox (43-58) amide can be used to translocate covalently attached peptides through biological membranes. The fluorescence label is useful for fluorescence-activated cell sorting (FACS) analysis and fluorescence microscope studies.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA