
Product Data Sheet

Product Name: Fmoc-4-cyano-Phe-OH

Cat. No.: GA21486

Chemical Properties

Cas. No. 173963-93-4

Formula $C_{25}H_{20}N_2O_4$ M.Wt 412.45

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Derivative for introducing the fluorophore p-cyanophenylalanine during Fmoc-SPPS. Educt for the synthesis of Fmoc-p-amidino-phenylalanine, which can be used for introducing an arginine mimetic.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA