
Product Data Sheet

Product Name: Fmoc-4-phenyl-D-Phe-OH

Cat. No.: GA21490

Chemical Properties

Cas. No. 205526-38-1

Formula $C_{30}H_{25}NO_4$

M.Wt 463.43

Solubility DMSO : 100 mg/mL (215.74 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA