

---

**Product Data Sheet**

---

Product Name: Fmoc-Arg(Me)<sub>2</sub>-OH (asymmetrical)

Cat. No.: GA21516

**Chemical Properties**

Cas. No. 268564-10-9

Formula C<sub>23</sub>H<sub>28</sub>N<sub>4</sub>O<sub>4</sub> M.Wt 424.5

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-ADMA, a building block for methylated histone sequences. The salt formation serves as temporary protection of the guanidino moiety.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA