
Product Data Sheet

Product Name: Fmoc-Arg(Me)₂-OH (symmetrical)

Cat. No.: GA21517

Chemical Properties

Cas. No. 823780-66-1

Formula C₂₃H₂₈N₄O₄ M.Wt 424.5

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-SDMA

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA