
Product Data Sheet

Product Name: Fmoc-Dap(Fmoc)-OH

Cat. No.: GA21587

Chemical Properties

Cas. No. 201473-90-7

Formula $C_{33}H_{28}N_2O_6$ M.Wt 548.6

Solubility DMSO : 100 mg/mL (182.29 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Suitable derivative for the Fmoc-SPPS of peptide dendrimers.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA