
Product Data Sheet

Product Name: Fmoc-D-His(Fmoc)-OH

Cat. No.: GA21621

Chemical Properties

Cas. No. 200926-18-7

Formula $C_{36}H_{29}N_3O_6$ M.Wt 599.64Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Mix of predominantly formed τ -Fmoc isomer and π -Fmoc isomer. CAS Number (τ -isomer) 200926-18-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA