
Product Data Sheet

Product Name: Fmoc-L-cysteic acid . disodium salt

Cat. No.: GA21706

Chemical Properties

Cas. No. 320384-09-6

Formula $C_{18}H_{15}NNa_2O_7S$ M.Wt 435.37

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Protected derivative for the introduction of this highly polar amino acid. Cya-containing peptides are obtained by oxidation of cysteine or cystine, as a post-translational modification, as a by-product during peptide synthesis, or, intentionally (by u

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA