
Product Data Sheet

Product Name: Fmoc-N-Me-Arg(Pbf)-OH

Cat. No.: GA21745

Chemical Properties

Cas. No. 913733-27-4

Formula C₃₅H₄₂N₄O₇S M.Wt 662.81

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Derivative for synthesizing peptides containing N α -methylated arginine, superior to Fmoc-N-Me-Arg(Mtr)-OH (4026642), as Pbf can be removed under much milder conditions.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA