
Product Data Sheet

Product Name: Fmoc-Orn(pyrazinylcarbonyl)-OH

Cat. No.: GA21761

Chemical Properties

Cas. No. 201046-61-9

Formula C₂₅H₂₄N₄O₅ M.Wt 460.49

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This ornithine derivative can be used as an effective substitute for the expensive cis-β-(4-pyrazinylcarbonylamino)cyclohexyl)alanine in the sequence of the two LHRH antagonists Argtide and Lystide. The antiovolatory activity of these antagonists r

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA