

---

**Product Data Sheet**

---

Product Name: Fmoc-Pro-Pro-Pro-Pro-Pro-OH

Cat. No.: GA21788

**Chemical Properties**

Cas. No. 454693-94-8

Formula C<sub>40</sub>H<sub>47</sub>N<sub>5</sub>O<sub>8</sub> M.Wt 725.84

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-pentaproline, Fmoc-PPPPP, potential building block for polyproline-based chiral stationary phases.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA