

---

**Product Data Sheet**

---

Product Name: Fmoc-Ser(PO<sub>3</sub>H<sub>2</sub>)-OH

Cat. No.: GA21795

**Chemical Properties**

Cas. No. 158171-15-4

Formula C<sub>18</sub>H<sub>18</sub>NO<sub>8</sub>P M.Wt 407.32

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Can be used for solid-phase synthesis of phosphopeptides by the Fmoc strategy.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA