
Product Data Sheet

Product Name: Fmoc-Tyr(SO₃H)-OH

Cat. No.: GA21831

Chemical Properties

Cas. No. 181952-24-9

Formula C₂₄H₂₁NO₈S M.Wt 483.5

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-sulfotyrosine can form the monosodium (CAS-No. 106864-37-3) or the disodium salt. A mix of both can be obtained as well. The salt form of Fmoc-Tyr(SO₃H)-OH varies from batch to batch.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA