
Product Data Sheet

Product Name: Galanin-Like Peptide (porcine)

Cat. No.: GA21891

Chemical Properties

Cas. No. 245114-96-9

Formula $C_{281}H_{443}N_{81}O_{78}$

M.Wt

6204.11

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

This member of the galanin family has been isolated from porcine hypothalamus. It is a 60 amino acid peptide that, unlike galanin, has a non-amidated C-terminus. GALP amino acid residues (9-21) are identical to the biologically active N-terminal (1-13)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA