
Product Data Sheet

Product Name: Galanin-Like Peptide (rat)

Cat. No.: GA21892

Chemical Properties

Cas. No. 245114-97-0

Formula C₂₈₈H₄₆₁N₈₇O₈₃S

M.Wt

6502.43

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Galanin-like peptide (GALP), first isolated from the porcine hypothalamus and subsequently also identified in rats and humans, is a 60 amino acid neuropeptide that unlike galanin, has a non-amidated C-terminus. Its amino acid residues (9-21) are identi

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA