
Product Data Sheet

Product Name: GHRP-2
Cat. No.: GA21901

Chemical Properties

Cas. No. 158861-67-7

Formula C₄₅H₅₅N₉O₆ M.Wt 817.99

Solubility >12.45 mg/mL in H₂O with ultrasonic; >74.1 mg/mL in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

GHRP-2 was studied in rats and produced a large plasma GH response in the presence of exogenous GH-releasing factor. This response was significantly augmented when somatostatin-antiserum was present. The hexapeptide was tested as diagnostic aid for GH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA