

---

**Product Data Sheet**

---

Product Name: Ginseng Tetrapeptide

Cat. No.: GA21902

**Chemical Properties**

Cas. No. 178553-95-2

Formula  $C_{18}H_{33}N_7O_7$  M.Wt 459.5

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Ginseng tetrapeptide stimulated the growth of baby hamster kidney (BHK)-21 cells, normal female bronchiogenic epithelioid (CCD-14 Br) cells, and normal human epidermal fibroblasts (NB1 RGB) without any morphological change and showed no toxicity.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA