
Product Data Sheet

Product Name: H-Aib-OBzl . HCl

Cat. No.: GA21959

Chemical Properties

Cas. No. 55456-40-1

Formula $C_{11}H_{15}NO_2 \cdot HCl$

M.Wt 229.71

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Building block for the solution synthesis of Aib-containing peptides, e.g. peptaibols. Coupling to this sterically demanding amino acid derivative allows to assess the efficiency of coupling reagents and protocols.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
