
Product Data Sheet

Product Name: H-Ala-Ala-Pro-pNA . HCl

Cat. No.: GA21980

Chemical Properties

Cas. No. 201732-27-6

Formula $C_{17}H_{23}N_5O_5 \cdot HCl$

M.Wt 413.86

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

AAP-pNA, chromogenic substrate for dipeptidyl aminopeptidase yscV, a thermosensitive enzyme from yeast and for prolyl tripeptidyl peptidases.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA