

---

**Product Data Sheet**

---

Product Name: H-Ala-His-Lys-OH

Cat. No.: GA22006

**Chemical Properties**

Cas. No. 126828-32-8

Formula  $C_{15}H_{26}N_6O_4$  M.Wt 354.41

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The tripeptide AHK (tripeptide-3), a GHK analog, also forms complexes with Cu(II). AHK-Cu promoted the growth of human hair follicles in vitro.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA