
Product Data Sheet

Product Name: H-Ala-Leu-Pro-Met-His-Ile-Arg-OH

Cat. No.: GA22011

Chemical Properties

Cas. No. 132160-04-4

Formula $C_{37}H_{64}N_{12}O_8S$ M.Wt 837.06Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The lactokinin ALPMHIR is an angiotensin-converting enzyme (ACE) inhibitory peptide, which corresponds to a sequence released by tryptic digestion from the milk protein β -lactoglobulin. It might be useful in the treatment of hypertension although

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
