
Product Data Sheet

Product Name: H-Ala-Phe-Lys-AMC

Cat. No.: GA22019

Chemical Properties

Cas. No. 120928-02-1

Formula C₂₈H₃₅N₅O₅ M.Wt 521.62

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

AFK-AMC, fluorogenic substrate for a simple and direct assay of plasmin, even in the presence of urokinase and thrombin.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA