

---

**Product Data Sheet**

---

Product Name: H-Arg(Me)-OH . HCl

Cat. No.: GA22047

**Chemical Properties**

Cas. No. 156706-47-7

Formula  $C_7H_{16}N_4O_2 \cdot HCl$ 

M.Wt 224.69

Solubility 25 mg/ml in methanol

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Cell-permeable inhibitor of all three isoforms of nitric oxide synthase (nNOS, eNOS, and iNOS). For the free base see H-Arg(Me)-OH.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA