
Product Data Sheet

Product Name: H-Arg(Me)-OH

Cat. No.: GA22048

Chemical Properties

Cas. No. 17035-90-4

Formula $C_7H_{16}N_4O_2$

M.Wt

188.23

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Cell-permeable inhibitor of all three isoforms of nitric oxide synthase: Ki 650 nM (nNOS), 700 nM (eNOS), 9500 nM (iNOS). The D-enantiomer can be used as negative control. CAS-No. (L-NMMA acetate): 53308-83-1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
