
Product Data Sheet

Product Name: H-Arg-Arg-Arg-Arg-Arg-Arg-OH

Cat. No.: GA22062

Chemical Properties

Cas. No. 96337-25-6

Formula C₃₆H₇₄N₂₄O₇ M.Wt 955.14

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure H-Arg-Arg-Arg-Arg-Arg-Arg-OH**Background**

Six arginines is the minimum number of residues required for obtaining a polycationic peptide which can effectively penetrate cells.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA