

---

**Product Data Sheet**

---

Product Name: H-Arg-Arg-Arg-Arg-OH

Cat. No.: GA22064

**Chemical Properties**

Cas. No. 26791-46-8

Formula C<sub>24</sub>H<sub>50</sub>N<sub>16</sub>O<sub>5</sub> M.Wt 642.77

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Tetraarginine, RRRR, is one of the commonly observed repeat units in protamine molecules. For the trifluoroacetate see H-Arg-Arg-Arg-Arg-OH.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA