
Product Data Sheet

Product Name: H-Arg-Phe-NH₂

Cat. No.: GA22089

Chemical Properties

Cas. No. 119051-99-9

Formula C₁₅H₂₄N₆O₂ M.Wt 320.4

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

RFamide has been used for raising antibodies employed in ELISA assays for detecting RFamide-type peptides in biological material.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA