
Product Data Sheet

Product Name: H-Asp(Phe-OH)-OH

Cat. No.: GA22120

Chemical Properties

Cas. No. 13433-10-8

Formula $C_{13}H_{16}N_2O_5$ M.Wt 280.28

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

β -Asp-Phe was isolated and identified in urine and plasma from healthy human volunteers. It has also been used in studies investigating structure-taste relationships of dipeptides.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA