
Product Data Sheet

Product Name: H-D-Ile-Phe-Lys-pNA

Cat. No.: GA22231

Chemical Properties

Cas. No. 76626-41-0

Formula C₂₇H₃₈N₆O₅ M.Wt 526.64

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Patented by Chromogenix AB, Sweden (CH 622 285). A highly specific substrate for human plasmin. For rat plasmin, Heinemann et al. determined a Km of 18.3 μM. Also cleaved by transmembrane tryptase (TMT).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA