
Product Data Sheet

Product Name: H-D-Phe-Phe-Arg-chloromethylketone

Cat. No.: GA22295

Chemical Properties

Cas. No. 74392-49-7

Formula $C_{25}H_{33}ClN_6O_3$ M.Wt 501.03

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Kallikrein inhibitor which is even effective in preventing apolipoprotein B-100 proteolysis in dextran sulfate (DS)-precipitated human plasma low density lipoproteins (LDL). D-FFR-CMK also inhibits gingipain R.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA