
Product Data Sheet

Product Name: H-D-Tyr-Val-Gly-OH

Cat. No.: GA22315

Chemical Properties

Cas. No. 86030-52-6

Formula $C_{16}H_{23}N_3O_5$ M.Wt 337.38

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The tripeptide H-D-Tyr-Val-Gly-OH can be converted into H-D-Tyr-Val-NH by a secretory granule-associated enzyme activity. Substrate (in combination with ascorbate) for peptidylglycine monooxygenase.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA