
Product Data Sheet

Product Name: H-D-Val-OBzl . p-tosylate

Cat. No.: GA22321

Chemical Properties

Cas. No. 17662-84-9

Formula $C_{12}H_{17}NO_2 \cdot C_7H_8O_3S$

M.Wt 379.48

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Useful building block for introduction of D-valine at the carboxy terminus in conventional peptide synthesis.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA