
Product Data Sheet

Product Name: Hepcidin-24 (human)

Cat. No.: GA22330

Chemical Properties

Cas. No.

Formula C₁₀₉H₁₆₅N₃₃O₂₈S₉

M.Wt 2674.31

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Hepcidin-24 is an excellent alternative to biotinylated or heavy isotope-labeled hepcidin-25 used as internal reference in the determination of Hep-25 by mass spectrometry-based methods, e.g. in human serum or urine.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA