
Product Data Sheet

Product Name: H-Ile-Arg-Pro-OH

Cat. No.: GA22562

Chemical Properties

Cas. No. 149756-66-1

Formula $C_{17}H_{32}N_6O_4$ M.Wt 384.48

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The tripeptide IRP showed ACE-inhibitory activity, IC50 1.8 μmol/l.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA