

---

**Product Data Sheet**

---

Product Name: H-Ile-Val-OH

Cat. No.: GA22581

**Chemical Properties**

Cas. No. 41017-96-3

Formula  $C_{11}H_{22}N_2O_3$  M.Wt 230.31

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

An effector dipeptide that enhances the catalytic activity of the proenzyme trypsinogen. Forms nanotubes as its retroanalog H-Val-Ile-OH. These pores absorb and store gases as methane or hydrogen.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA