
Product Data Sheet

Product Name: Hippuryl-Phe-Arg-OH

Cat. No.: GA22588

Chemical Properties

Cas. No. 73167-83-6

Formula $C_{24}H_{30}N_6O_5$ M.Wt 482.54

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Hip-Phe-Arg is used as substrate for angiotensin I-converting enzyme and for dipeptidyl carboxypeptidase.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA