
Product Data Sheet

Product Name: Histone H3 (1-21)-Gly-Gly-Lys(biotinyl) amide

Cat. No.: GA22595

Chemical Properties

Cas. No. 2243207-10-3

Formula $C_{114}H_{205}N_{43}O_{32}S$

M.Wt 2722.22

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA